






















CHOOSE YOUR SET

Choose the all-in-one set that best suits your needs, according to your veterinary practice profile and care level, or according to the treatments you usually practice.

		SET LEVEL 1	SET LEVEL 2	SET LEVEL 3	SET LEVEL 4
			 	 	 
Motor & Drive Shaft	 Motor 1 Drive shaft 47 inches	✓			
	 Motor EVOLUTION - Drive shaft 47 inches		✓	✓	✓
	 Flexidisc curved long 21.3 inches	✓	✓	✓	✓
	 06 profiling handpiece		✓	✓	✓
Handpiece	 Polyfloat straight long Watercooling 21.3 inches			✓	✓
	 Polyfloat applecore diam. 14			✓	✓
	 Polyfloat occlusal burr				✓
	 Polyfloat conical burr				✓
	 Polyfloat long versatile burr				✓
	 11 th hand piece				✓
	 Stainless Steel Brush with support	✓	✓	✓	✓
Instrument Case	 Instrument case	✓			
	 Trolley		✓	✓	✓

	SET LEVEL 1	SET LEVEL 2	SET LEVEL 3	SET LEVEL 4
 RETAIL PRICE * with watercooling**	not available	\$ 7,500.00	\$ 10,000.00	\$ 11,500.00
		J1410s2wc	J1410s3wc	J1410s4wc
RETAIL PRICE * without watercooling	\$ 3,500.00	\$ 6,500.00	\$ 8,500.00	\$ 10,000.00
	J1411s1	J1410s2	J1410s3	J1410s4

* Retail price excluding transport costs and taxes valid until 03.31.2020 - ** Watercooling is available exclusively for Flexidisc and Polyfloat Handpieces

Flexxidisc Straight Short

[J1413d1]



- Perfect for incisors and 06th profiling
- Length: 7.09 inches (18 cm)
- Disc diameter: 1.26 inch (32 mm)

Polyfloat Burrs

[J1412d9]



Ø 0.55 inch (14 mm)

[J1412d8]



Ø 0.70 inch (18 mm)

[J1412d11]



Ø 0.70 inch (18 mm)

- Perfect for incisors, canine and 06th profiling

- Recommended for floating finishing

Ergonomic Handle

[J1414d10 - J1414d7a]



- Improves grip for short handpieces (06th profiling)
- Very easy to install

Lithium Ion Battery

[J1410d3]



- Autonomy: up to 10 horses
- Complete loading time: 2h30
- Weight: 815 g

Evolution Battery Car Charger

[J1410d5]



- Perfect for charging the battery while travelling

Right or Left Holster

[J1410hr - J1410hL]



- Safely keeps the equipment closeby
- Eases the movements

HDE NORTH AMERICA LLC

11999 Katy Freeway, #330, Houston TX 77079
Cell: (713) 444-2531 - contact@horse-dental-equipment.us

www.horse-dental-equipment.us