

HYDRO JET II

The Hydrojet II is the essential tool for every equine dental practitioner interested in a versatile instrument for the efficient flushing of diastemas as well as the rinsing of the mouth.



The Hydro jet II includes the following:

- Hydro Jet II with battery and charger
- Rinsing nozzle
- Diastema nozzle
- 29.5-ft / 6-meter water tube
- 0.26 gal US / 1-liter bottle
- Slim carriage bag

ACC_SHJ2

From US \$771

The Hydro jet II has been completely redesigned for easier work with several new smart features:



Fast & easy brass nozzle connectors



Water flow regulation from 0.22 to 0.45 gal US/min - 0,75 to 1,5 liter/minute



Battery level indicator



Curved rinsing nozzle



Rounded-off diastema tip



Battery moved on the top

- Waterpressure: 7 bars
- Battery usage time: 2 hours of continuous use
- Battery charging time: 1,5 - 2 hours
- The diastema nozzle is perfect for flushing diastemas and removing interdental food residues. Its length and shape allow easy access in the back of the mouth.
- The rinsing nozzle is for general mouth cleaning all along the different dental treatment stages. Its specific curved shape naturally fits the mouth.