



HORSE DENTAL  
EQUIPMENT

# EQUUS JET RANGE



# EQUUS JET RANGE

## Multi-purpose rinsing devices

The Equus Jet range has been developed as an autonomous, robust, compact cleaning and cooling system with the mobile equine dentistry practitioners and clinics in mind.

It is ideal for cleaning the mouth before and after dental work, diastema, rinsing after tooth extraction, sinus flushing or cooling when working with handpieces.

## EQUUS JET ECO

Product features:

- High pressure (up to 10 bars)
- Flow regulated via the wands
- Lead battery (up to a week autonomy)
- Fast connecting and disconnecting tubes and wands
- Double filter system (external and internal)



The Equus Jet Eco includes the following accessories:

- One flushing wand
- One diastema wand
- An input polyurethane tube with a strainer and filter
- An output spiral polyurethane tube with fast connector on one side and fast connector with valve on the other, allowing to change the wands without shutting off the pump to save time
- A special lead battery charger

Art.-No. #	Description
PEGA_EJ-EJE02	Equus Jet Eco
PEGA_EJ-CBEP01	Equus Jet Eco charger
PEGA_EJ-LDEJ1	Diastema handpiece
PEGA_EJ-LDEJ2	Rinsing handpiece

# EQUUS JET

Product features:

- High pressure (up to 10 bars)
- 5 set levels of pressure regulated via the white control panel on top
- Lithium ion battery (up to 3 weeks autonomy)
- Possibility to use a wired pedal to gradually regulate flow and pressure
- Fast connecting and disconnecting tubes and wands
- Double filter system (external and internal)



The Equus Jet includes the following accessories:

- One flushing wand
- One diastema wand
- An input polyurethane tube with strainer and filter
- An output spiral polyurethane tube with fast connector on one side and fast connector with valve on the other, allowing to change the wands without shutting off the pump to save time
- A special lithium ion battery charger

Art.-No. #	Description
PEGA_EJH-EJ11	Equus Jet
PEGA_EJ-EJ04	Equus Jet - Option Pegasos Trolley option
PEGA_EJ-CBLI14	Equus Jet charger
PEGA_EJ-LDEJ1	Diastema handpiece
PEGA_EJ-LDEJ2	Rinsing handpiece
PEGA_EJH-PDC02	Equus Jet Foot Pedal
PEGA_EJH-KLS02	Equus Jet Sinus Kit

# EQUUS JET CONNECT

## Product features:

- High pressure (up to 10 bars)
- Control screen with general settings and display : pressure, temperature, battery level, pumps, pedal, pairing
- Bluetooth and Android application
- Lithium ion battery (up to 3 weeks autonomy)
- Possibility to use a wired pedal to gradually regulate flow and pressure
- Fast connecting and disconnecting tubes and wands
- Double filter system (external and internal)



## The Equus Jet Connect the following accessories:

- One flushing wand
- One diastema wand
- An input polyurethane tube with a strainer and filter
- An output spiral polyurethane tube with fast connector on one side and fast connector with valve on the other, allowing to change the wands without shutting off the pump to save time
- An outlet tube with fast connector on the pump side for the cooling pump
- An tube with luer connector for the peristaltic pump
- A special lithium ion battery charger

Art.-No. #	Description
PEGA_EJ-EJV21-2	Equus Jet Connect with 2 pumps (high pressure rinsing and low pressure watercooling)
PEGA_EJ-CBLI14	Equus Jet charger
PEGA_EJ50630	Peristaltic pump
PEGA_EJH-PDC02	Equus Jet Foot Pedal
PEGA_EL-PC01	Wireless Pedal
PEGA_EJH-KLS02	Equus Jet Sinus Kit
PEGA_EJ-LDLL90	Diastema luer lock wand for Equus Jet
PEGA_EJ-AIDG18	Luer lock Needles Box of 10 needles in several lengths and diameters
PEGA_EJ-EDB02	Dimmer Box Connect

# EQUUS JET RANGE

## Comparative table

	EQUUS JET ECO	EQUUS JET	EQUUS JET CONNECT
Battery type	Lead	Lithium ion	Lithium ion
Battery autonomy**	Up to a week**	Up to 3 weeks**	Up to 3 weeks**
Pressure level	Up to 10 bars	Up to 10 bars	Up to 10 bars
Number of pressure levels	1	5	6
Control panel with pressure regulation		✓	✓
Fast connecting and disconnecting tubes and wands	✓	✓	✓
Double filter system (external and internal)	✓	✓	✓
High pressure pump	✓	✓	✓
Cooling pump			✓
Peristaltic pump			✓ *
Pegasos trolley option		✓ *	✓ *
Control screen with general settings and display : pressure, temperature, battery level, pumps, pedal, pairing			✓
Bluetooth and Android application			✓
Wired foot pedal control		✓ *	
Wireless foot pedal control			
Kit sinus		✓ *	✓ *
			

\* Optional. \*\*Average autonomy.

# EQUUS JET RANGE ACCESSORIES

## EQUUS JET ECO CHARGER

Different sockets available:  
French, Swiss, UK, USA



## EQUUS JET AND EQUUS JET CONNECT CHARGER

Different sockets available:  
French, Swiss, UK, USA



## WIRED FOOT PEDAL

The wired foot pedal allows you to control the pressure of the Equus Jet directly with your foot.

You keep freedom of movement as well as the availability of your hands while keeping the pressure setting.

This pedal is used with Equus Jet.



## WIRELESS FOOT PEDAL CONNECT

This Bluetooth wireless Foot Pedal Connect, allows you to control:

- the speed of the Equus Jet Connect
- the Equus Jet Connect pump
- the Dimmer Box Connect
- the options such as the reversal the motor rotation direction, the start or stop of the water pump running with the motor.



## PEGASOS TROLLEY OPTION

Thanks to this system you can fix your Equus Jet on your Pegasus trolley.



### RINSING WAND

Rinse the horse's mouth in one single gesture, with upper and lower opening. High flow with slight pressure to remove food debris before and after equine dental treatment. Its length allows you to easily access the back of the mouth without hurting the horse.



### DIASTEMA WAND

Diastema wand with a 90° angle to clean interdental spaces, diastemas and cavities after extraction, wounds. Its length allows you to easily access the back of the mouth without hurting the horse.



### DIASTEMA LUER LOCK WAND

This new diastema lance allows you to use needles of different lengths as well as to tilt them in the desired area.



### LUER LOCK NEEDLES

Box of 10 straight stainless needles for Equus Jet luer lock wand in several lengths and diameters.



### SINUS KIT

This kit allows you to clean the horse's sinuses with the Equus Jet. It's recommended to use the pedal to vary the pressure flow.



### PERISTALTIC PUMP



### DIMMER BOX CONNECT

Equus Dimmer Box Connect, connected management box for the use of a wireless foot pedal and to control the Storz endoscope.

